

# Film Array Standard Standard for Polariscope

It is available for color sample of retardation, and also, for standard piece to calibrate a Polariscope.

## FAS200

List price JPY100,000



- ❖ Convenient to quantify retardation in a sample
- ❖ Available to watch color gradation by retardation on Sensitive tint color method
- ❖ Including with the certification of test.

### Specification

Measurement range	: 0-200nm
Accuracy	: ±5nm
Material	: Film
Size	: 45x210x10mm
Weight	: 100g

The retardation values on above photo are reference. Each product has own values.

Orihara Industrial Co., Ltd.  
5-47-15, Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013, Japan  
TEL +81-3-3985-9531 FAX +81-3-3985-9532  
<http://www.orihara-ss.co.jp> e-mail [orihara@orihara-ss.co.jp](mailto:orihara@orihara-ss.co.jp)